

# WATER TEST

## A New Method For The Detection Of Coliforms-enzyme Substrate Method

Coliform group / Fecal coliforms / heat-resistant coliforms / Escherichia coli

### METHOD ADVANTAGE

1. No need to operate in a sterile room.
2. The manual operation time is less than 1 minute.
3. No need for medium preparation and bulk glassware sterilization.
4. Qualitative and quantitative analysis can be completed within 24 hours without verification test.




**LIFECOSM BIOTECH LIMITED**



# THE PRINCIPLE OF IMMOBILIZED SUBSTRATE TECHNOLOGY ENZYME SUBSTRATE METHOD

The principle of the enzyme substrate method uses the total coliform to produce  $\beta$ -galactosidase to decompose the color source substrate-ONPG (Ortho-nitrophenyl- $\beta$ -D-galactopyranoside) to make the culture medium yellow; Escherichia coli produces  $\beta$ -glucuronidase. The principle that aldolase ( $\beta$ -glucuronidase) decomposes the fluorescent color source substrate MUG (4-methyl-umbelliferyl- $\beta$ -glucuronide) to make the culture solution generate fluorescence under ultraviolet light with a wavelength of 366 nm, to qualitatively quantify the total coliform bacteria in water, Fecal coliforms (heat-resistant coliforms) and Escherichia coli.

Lifecosc enzyme substrate method reagent is Minimal Medium ONPG-MUG (MMO-MUG) and HJ 1001- 2018's medium. When the Lifecosc enzyme substrate method detection reagent is added to the water sample containing the total coliform bacteria group, the target bacteria are propagated in the medium containing the immobilized substrate Minimal Medium ONPG-MUG, and the specific biological enzyme  $\beta$  produced by the coliform group bacteria group -Galactosidase can decompose the fixed chromogen substrate ONPG in Minimal Medium ONPG-MUG medium, making the medium yellow; at the same time, Escherichia coli produces specific  $\beta$ -glucuronidase in water samples to decompose Minimal Medium The immobilized fluorescent substrate MUG in ONPG-MUG medium produces characteristic blue fluorescence. The same principle, heat-resistant coliform (fecal coliform) will decompose the chromogenic substrate ONPG in Minimal Medium ONPG-MUG medium when cultured at 44.5°C, making the medium appear yellow. Therefore, it is possible to qualitatively detect total coliforms, heat-resistant (fecal) coliforms and Escherichia coli in water. When the biological enzymes of a few non-target bacteria in the water sample try to react with ONPG or MUG and interfere with the test results, the specially formulated matrix in the Lifecosc medium will inhibit the reaction process and control the generation of false positives; at the same time, these specially formulated matrices can It effectively promotes the reaction between the biological enzymes of the target bacteria and the immobilized substrate, so the incidence of false negatives is lower than that of traditional methods. Lifecosc enzyme substrate method reagent is used in combination with 51 or 97-well quantitative detection disk to quantitatively detect target bacteria in 100 ml of water samples. without dilution.



# ENZYME SUBSTRATE ASSAY

## QUALITATIVE DETECTION



1

Add reagent to 100ml water sample, after dissolving, incubate at 36°C for 24h



2

Interpretation of results:

colorless = negative

Yellow = positive for total coliforms

Yellow + fluorescence = Escherichia coli positive.

## QUANTITATIVE DETECTION



1

Add the reagents to the water sample and mix well.



2

Pour into 51-well quantitative detection plate (quantitative well plate) or 97-well quantitative detection plate (quantitative well plate)



3

Use the program-controlled quantitative sealing machine to seal the quantitative detection disc (quantitative well plate) for sealing and incubate at 36°C for 24h

Heat-resistant coliform/fecal coliform culture at 44.5°C for 24h is yellow and positive



4

Interpretation of results:

colorless = negative

Yellow checkered = positive total coliforms

Yellow + fluorescent grid = Escherichia coli positive  
reference MPN table count



**97-WELL QUANTITATIVE PLATE MPN TABLE (100)**

Large Wells Positive	SMALL WELLS POSITIVE																								
	0	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
0	<1	1.0	2.0	3.0	4.0	5.0	6.0	7.0	8.0	9.0	10.0	11.0	12.0	13.0	14.1	15.1	16.1	17.1	18.1	19.1	20.2	21.2	22.2	23.3	24.3
1	1.0	2.0	3.0	4.0	5.0	6.0	7.1	8.1	9.1	10.1	11.1	12.1	13.2	14.2	15.2	16.2	17.3	18.3	19.3	20.4	21.4	22.4	23.5	24.5	25.6
2	2.0	3.0	4.1	5.1	6.1	7.1	8.1	9.2	10.2	11.2	12.2	13.3	14.3	15.4	16.4	17.4	18.5	19.5	20.6	21.6	22.7	23.7	24.8	25.8	26.9
3	3.1	4.1	5.1	6.1	7.2	8.2	9.2	10.3	11.3	12.4	13.4	14.5	15.5	16.5	17.6	18.6	19.7	20.8	21.8	22.9	23.9	25.0	26.1	27.1	28.2
4	4.1	5.2	6.2	7.2	8.3	9.3	10.4	11.4	12.5	13.5	14.6	15.6	16.7	17.8	18.8	19.9	21.0	22.0	23.1	24.2	25.3	26.3	27.4	28.5	29.6
5	5.2	6.3	7.3	8.4	9.4	10.5	11.5	12.6	13.7	14.7	15.8	16.9	17.9	19.0	20.1	21.2	22.2	23.3	24.4	25.5	26.6	27.7	28.8	29.9	31.0
6	6.3	7.4	8.4	9.5	10.6	11.6	12.7	13.8	14.9	16.0	17.0	18.1	19.2	20.3	21.4	22.5	23.6	24.7	25.8	26.9	28.0	29.1	30.2	31.3	32.4
7	7.5	8.5	9.6	10.7	11.8	12.8	13.9	15.0	16.1	17.2	18.3	19.4	20.5	21.6	22.7	23.8	24.9	26.0	27.1	28.3	29.4	30.5	31.6	32.8	33.9
8	8.6	9.7	10.8	11.9	13.0	14.1	15.2	16.3	17.4	18.5	19.6	20.7	21.8	22.9	24.1	25.2	26.3	27.4	28.6	29.7	30.8	32.0	33.1	34.3	35.4
9	9.8	10.9	12.0	13.1	14.2	15.3	16.4	17.6	18.7	19.8	20.9	22.0	23.2	24.3	25.4	26.6	27.7	28.9	30.0	31.2	32.3	33.5	34.6	35.8	37.0
10	11.0	12.1	13.2	14.4	15.5	16.6	17.7	18.9	20.0	21.1	22.3	23.4	24.6	25.7	26.9	28.0	29.2	30.3	31.5	32.7	33.8	35.0	36.2	37.4	38.6
11	12.2	13.4	14.5	15.6	16.8	17.9	19.1	20.2	21.4	22.5	23.7	24.8	26.0	27.2	28.3	29.5	30.7	31.9	33.0	34.2	35.4	36.6	37.8	39.0	40.2
12	13.5	14.6	15.8	16.9	18.1	19.3	20.4	21.6	22.8	23.9	25.1	26.3	27.5	28.6	29.8	31.0	32.2	33.4	34.6	35.8	37.0	38.2	39.5	40.7	41.9
13	14.8	16.0	17.1	18.3	19.5	20.6	21.8	23.0	24.2	25.4	26.6	27.8	29.0	30.2	31.4	32.6	33.8	35.0	36.2	37.5	38.7	39.9	41.2	42.4	43.6
14	16.1	17.3	18.5	19.7	20.9	22.1	23.3	24.5	25.7	26.9	28.1	29.3	30.5	31.7	33.0	34.2	35.4	36.7	37.9	39.1	40.4	41.6	42.9	44.2	45.4
15	17.5	18.7	19.9	21.1	22.3	23.5	24.7	25.9	27.2	28.4	29.6	30.9	32.1	33.3	34.6	35.8	37.1	38.4	39.6	40.9	42.2	43.4	44.7	46.0	47.3
16	18.9	20.1	21.3	22.6	23.8	25.0	26.2	27.5	28.7	30.0	31.2	32.5	33.7	35.0	36.3	37.5	38.8	40.1	41.4	42.7	44.0	45.3	46.6	47.9	49.2
17	20.3	21.6	22.8	24.1	25.3	26.6	27.8	29.1	30.3	31.6	32.9	34.1	35.4	36.7	38.0	39.3	40.6	41.9	43.2	44.5	45.9	47.2	48.5	49.8	51.2
18	21.8	23.1	24.3	25.6	26.9	28.1	29.4	30.7	32.0	33.3	34.6	35.9	37.2	38.5	39.8	41.1	42.4	43.8	45.1	46.5	47.8	49.2	50.5	51.9	53.2
19	23.3	24.6	25.9	27.2	28.5	29.8	31.1	32.4	33.7	35.0	36.3	37.6	39.0	40.3	41.6	43.0	44.3	45.7	47.1	48.4	49.8	51.2	52.6	54.0	55.4
20	24.9	26.2	27.5	28.8	30.1	31.5	32.8	34.1	35.4	36.8	38.1	39.5	40.8	42.2	43.6	44.9	46.3	47.7	49.1	50.5	51.9	53.3	54.7	56.1	57.6
21	26.5	27.9	29.2	30.5	31.8	33.2	34.5	35.9	37.3	38.6	40.0	41.4	42.8	44.1	45.5	46.9	48.4	49.8	51.2	52.6	54.1	55.5	56.9	58.4	59.9
22	28.2	29.5	30.9	32.3	33.6	35.0	36.4	37.7	39.1	40.5	41.9	43.3	44.8	46.2	47.6	49.0	50.5	51.9	53.4	54.8	56.3	57.8	59.3	60.8	62.3
23	29.9	31.3	32.7	34.1	35.5	36.8	38.3	39.7	41.1	42.5	43.9	45.4	46.8	48.3	49.7	51.2	52.7	54.2	55.6	57.1	58.6	60.2	61.7	63.2	64.7
24	31.7	33.1	34.5	35.9	37.3	38.8	40.2	41.7	43.1	44.6	46.0	47.5	49.0	50.5	52.0	53.5	55.0	56.5	58.0	59.5	61.1	62.6	64.2	65.8	67.3
25	33.6	35.0	36.4	37.9	39.3	40.8	42.2	43.7	45.2	46.7	48.2	49.7	51.2	52.7	54.3	55.8	57.3	58.9	60.5	62.0	63.6	65.2	66.8	68.4	70.0
26	35.5	36.9	38.4	39.9	41.4	42.8	44.3	45.9	47.4	48.9	50.4	52.0	53.5	55.1	56.7	58.2	59.8	61.4	63.0	64.7	66.3	67.9	69.6	71.2	72.9
27	37.4	38.9	40.4	42.0	43.5	45.0	46.5	48.1	49.6	51.2	52.8	54.4	56.0	57.6	59.2	60.8	62.4	64.1	65.7	67.4	69.1	70.8	72.5	74.2	75.9
28	39.5	41.0	42.6	44.1	45.7	47.3	48.8	50.4	52.0	53.6	55.2	56.9	58.5	60.2	61.8	63.5	65.2	66.9	68.6	70.3	72.0	73.7	75.5	77.3	79.0
29	41.7	43.2	44.8	46.4	48.0	49.6	51.2	52.8	54.5	56.1	57.8	59.5	61.2	62.9	64.6	66.3	68.0	69.8	71.5	73.3	75.1	76.9	78.7	80.5	82.4
30	43.9	45.5	47.1	48.7	50.4	52.0	53.7	55.4	57.1	58.8	60.5	62.2	64.0	65.7	67.5	69.3	71.0	72.9	74.7	76.5	78.3	80.2	82.1	84.0	85.9
31	46.2	47.9	49.5	51.2	52.9	54.6	56.3	58.1	59.8	61.6	63.3	65.1	66.9	68.7	70.5	72.4	74.2	76.1	78.0	79.9	81.8	83.7	85.7	87.6	89.6
32	48.7	50.4	52.1	53.8	55.6	57.3	59.1	60.9	62.7	64.5	66.3	68.2	70.0	71.9	73.8	75.7	77.6	79.5	81.5	83.5	85.4	87.5	89.5	91.5	93.6
33	51.2	53.0	54.8	56.5	58.3	60.2	62.0	63.8	65.7	67.6	69.5	71.4	73.3	75.2	77.2	79.2	81.2	83.2	85.2	87.3	89.3	91.4	93.6	95.7	97.8
34	53.9	55.7	57.6	59.4	61.3	63.1	65.0	67.0	68.9	70.8	72.8	74.8	76.8	78.8	80.8	82.9	85.0	87.1	89.2	91.4	93.5	95.7	97.9	100.2	102.4
35	56.8	58.6	60.5	62.4	64.4	66.3	68.3	70.3	72.3	74.3	76.3	78.4	80.5	82.6	84.7	86.9	89.1	91.3	93.5	95.7	98.0	100.3	102.6	105.0	107.3
36	59.8	61.7	63.7	65.7	67.7	69.7	71.7	73.8	75.9	78.0	80.1	82.3	84.5	86.7	88.9	91.2	93.5	95.8	98.1	100.5	102.9	105.3	107.7	110.2	112.7
37	62.9	65.0	67.0	69.1	71.2	73.3	75.4	77.6	79.8	82.0	84.2	86.5	88.8	91.1	93.4	95.8	98.2	100.6	103.1	105.6	108.1	110.7	113.3	115.9	118.6
38	66.3	68.4	70.6	72.7	74.9	77.1	79.4	81.6	83.9	86.2	88.6	91.0	93.4	95.8	98.3	100.8	103.4	105.9	108.6	111.2	113.9	116.6	119.4	122.2	125.0
39	70.0	72.2	74.4	76.7	78.9	81.3	83.6	86.0	88.4	90.9	93.4	95.9	98.4	101.0	103.6	106.3	109.0	111.8	114.6	117.4	120.3	123.2	126.1	129.2	132.2
40	73.8	76.2	78.5	80.9	83.3	85.7	88.2	90.8	93.3	95.9	98.5	101.2	103.9	106.7	109.5	112.4	115.3	118.2	121.2	124.3	127.4	130.5	133.7	137.0	140.3
41	78.0	80.5	83.0	85.5	88.0	90.6	93.3	95.9	98.7	101.4	104.3	107.1	110.0	113.0	116.0	119.1	122.2	125.4	128.7	132.0	135.4	138.8	142.3	145.9	149.5
42	82.6	85.2	87.8	90.5	93.2	96.0	98.8	101.7	104.6	107.6	110.6	113.7	116.9	120.1	123.4	126.7	130.1	133.6	137.2	140.8	144.5	148.3	152.2	156.1	160.2
43	87.6	90.4	93.2	96.0	99.0	101.9	105.0	108.1	111.2	114.5	117.8	121.1	124.6	128.1	131.7	135.4	139.1	143.0	147.0	151.0	155.2	159.4	163.8	168.2	172.8
44	93.1	96.1	99.1	102.2	105.4	108.6	111.9	115.3	118.7	122.3	125.9	129.6	133.4	137.4	141.4	145.5	149.7	154.1	158.5	163.1	167.9	172.7	177.7	182.9	188.2
45	99.3	102.5	105.8	109.2	112.6	116.2	119.8	123.6	127.4	131.4	135.4	139.6	143.9	148.3	152.9	157.6	162.4	167.4	172.6	178.0	183.5	189.2	195.1	201.2	207.5
46	106.3	109.8	113.4	117.2	121.0	125.0	129.1	133.3	137.6	142.1	146.7	151.5	156.5	161.6	167.0	172.5	178.2	184.2	190.4	196.8	203.5	210.5	217.8	225.4	233.3
47	114.3	118.3	122.4	126.6	130.9	135.4	140.1	145.0	150.0	155.3	160.7	166.4	172.3	178.5	185.0	191.8	198.9	206.4	214.2	222.4	231.0	240.0	249.5	259.5	270.0
48	123.9	128.4	133.1	137.9	143.0	148.3	153.9	159.7	165.8	172.2	178.9	186.0	193.5	201.4	209.8	218.7	228.2	238.2	248.9	260.3	272.3	285.1	298.7	313.0	328.2
49	135.5	140.8	146.4	152.3	158.5	165.0	172.0	179.3	187.2	195.6	204.6	214.3	224.7	235.9	248.1	261.3	275.5	290.9	307.6	325.5					

**97-WELL QUANTITATIVE PLATE MPN TABLE (100)**




Large Wells Positive	SMALL WELLS POSITIVE																							
	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48
0	25.3	26.4	27.4	28.4	29.5	30.5	31.5	32.6	33.6	34.7	35.7	36.8	37.8	38.9	40.0	41.0	42.1	43.1	44.2	45.3	46.3	47.4	48.5	49.5
1	26.6	27.7	28.7	29.8	30.8	31.9	32.9	34.0	35.0	36.1	37.2	38.2	39.3	40.4	41.4	42.5	43.6	44.7	45.7	46.8	47.9	49.0	50.1	51.2
2	27.9	29.0	30.0	31.1	32.2	33.2	34.3	35.4	36.5	37.5	38.6	39.7	40.8	41.9	43.0	44.0	45.1	46.2	47.3	48.4	49.5	50.6	51.7	52.8
3	29.3	30.4	31.4	32.5	33.6	34.7	35.8	36.8	37.9	39.0	40.1	41.2	42.3	43.4	44.5	45.6	46.7	47.8	48.9	50.0	51.2	52.3	53.4	54.5
4	30.7	31.8	32.8	33.9	35.0	36.1	37.2	38.3	39.4	40.5	41.6	42.8	43.9	45.0	46.1	47.2	48.3	49.5	50.6	51.7	52.9	54.0	55.1	56.3
5	32.1	33.2	34.3	35.4	36.5	37.6	38.7	39.9	41.0	42.1	43.2	44.4	45.5	46.6	47.7	48.9	50.0	51.2	52.3	53.5	54.6	55.8	56.9	58.1
6	33.5	34.7	35.8	36.9	38.0	39.2	40.3	41.4	42.6	43.7	44.8	46.0	47.1	48.3	49.4	50.6	51.7	52.9	54.1	55.2	56.4	57.6	58.7	59.9
7	35.0	36.2	37.3	38.4	39.6	40.7	41.9	43.0	44.2	45.3	46.5	47.7	48.8	50.0	51.2	52.3	53.5	54.7	55.9	57.1	58.3	59.4	60.6	61.8
8	36.6	37.7	38.9	40.0	41.2	42.3	43.5	44.7	45.9	47.0	48.2	49.4	50.6	51.8	53.0	54.1	55.3	56.5	57.7	59.0	60.2	61.4	62.6	63.8
9	38.1	39.3	40.5	41.6	42.8	44.0	45.2	46.4	47.6	48.8	50.0	51.2	52.4	53.6	54.8	56.0	57.2	58.4	59.7	60.9	62.1	63.4	64.6	65.8
10	39.7	40.9	42.1	43.3	44.5	45.7	46.9	48.1	49.3	50.6	51.8	53.0	54.2	55.5	56.7	57.9	59.2	60.4	61.7	62.9	64.2	65.4	66.7	67.9
11	41.4	42.6	43.8	45.0	46.3	47.5	48.7	49.9	51.2	52.4	53.7	54.9	56.1	57.4	58.6	59.9	61.2	62.4	63.7	65.0	66.3	67.5	68.8	70.1
12	43.1	44.3	45.6	46.8	48.1	49.3	50.6	51.8	53.1	54.3	55.6	56.8	58.1	59.4	60.7	62.0	63.2	64.5	65.8	67.1	68.4	69.7	71.0	72.4
13	44.9	46.1	47.4	48.6	49.9	51.2	52.5	53.7	55.0	56.3	57.6	58.9	60.2	61.5	62.8	64.1	65.4	66.7	68.0	69.3	70.7	72.0	73.3	74.7
14	46.7	48.0	49.3	50.5	51.8	53.1	54.4	55.7	57.0	58.3	59.6	60.9	62.3	63.6	64.9	66.3	67.6	68.9	70.3	71.6	73.0	74.4	75.7	77.1
15	48.6	49.9	51.2	52.5	53.8	55.1	56.4	57.8	59.1	60.4	61.8	63.1	64.5	65.8	67.2	68.5	69.9	71.3	72.6	74.0	75.4	76.8	78.2	79.6
16	50.5	51.8	53.2	54.5	55.8	57.2	58.5	59.9	61.2	62.6	64.0	65.3	66.7	68.1	69.5	70.9	72.3	73.7	75.1	76.5	77.9	79.3	80.8	82.2
17	52.5	53.9	55.2	56.6	58.0	59.3	60.7	62.1	63.5	64.9	66.3	67.7	69.1	70.5	71.9	73.3	74.8	76.2	77.6	79.1	80.5	82.0	83.5	84.9
18	54.6	56.0	57.4	58.8	60.2	61.6	63.0	64.4	65.8	67.2	68.6	70.1	71.5	73.0	74.4	75.9	77.3	78.8	80.3	81.8	83.3	84.8	86.3	87.8
19	56.8	58.2	59.6	61.0	62.4	63.9	65.3	66.8	68.2	69.7	71.1	72.6	74.1	75.5	77.0	78.5	80.0	81.5	83.1	84.6	86.1	87.6	89.2	90.7
20	59.0	60.4	61.9	63.3	64.8	66.3	67.7	69.2	70.7	72.2	73.7	75.2	76.7	78.2	79.8	81.3	82.8	84.4	85.9	87.5	89.1	90.7	92.2	93.8
21	61.3	62.8	64.3	65.8	67.3	68.8	70.3	71.8	73.3	74.9	76.4	77.9	79.5	81.1	82.6	84.2	85.8	87.4	89.0	90.6	92.2	93.8	95.4	97.1
22	63.8	65.3	66.8	68.3	69.8	71.4	72.9	74.5	76.1	77.6	79.2	80.8	82.4	84.0	85.6	87.2	88.9	90.5	92.1	93.8	95.5	97.1	98.8	100.5
23	66.3	67.8	69.4	71.0	72.5	74.1	75.7	77.3	78.9	80.5	82.2	83.8	85.4	87.1	88.7	90.4	92.1	93.8	95.5	97.2	98.9	100.6	102.4	104.1
24	68.9	70.5	72.1	73.7	75.3	77.0	78.6	80.3	81.9	83.6	85.2	86.9	88.6	90.3	92.0	93.8	95.5	97.2	99.0	100.7	102.5	104.3	106.1	107.9
25	71.7	73.3	75.0	76.6	78.3	80.0	81.7	83.3	85.1	86.8	88.5	90.2	92.0	93.7	95.5	97.3	99.1	100.9	102.7	104.5	106.3	108.2	110.0	111.9
26	74.6	76.3	78.0	79.7	81.4	83.1	84.8	86.6	88.4	90.1	91.9	93.7	95.5	97.3	99.2	101.0	102.9	104.7	106.6	108.5	110.4	112.3	114.2	116.2
27	77.6	79.4	81.1	82.9	84.6	86.4	88.2	90.0	91.9	93.7	95.5	97.4	99.3	101.2	103.1	105.0	106.9	108.8	110.8	112.7	114.7	116.7	118.7	120.7
28	80.8	82.6	84.4	86.3	88.1	89.9	91.8	93.7	95.6	97.5	99.4	101.3	103.3	105.2	107.2	109.2	111.2	113.2	115.2	117.3	119.3	121.4	123.5	125.6
29	84.2	86.1	87.9	89.8	91.7	93.7	95.6	97.5	99.5	101.5	103.5	105.5	107.5	109.5	111.6	113.7	115.7	117.8	120.0	122.1	124.2	126.4	128.6	130.8
30	87.8	89.7	91.7	93.6	95.6	97.6	99.6	101.6	103.7	105.7	107.8	109.9	112.0	114.2	116.3	118.5	120.6	122.8	125.1	127.3	129.5	131.8	134.1	136.4
31	91.6	93.6	95.6	97.7	99.7	101.8	103.9	106.0	108.2	110.3	112.5	114.7	116.9	119.1	121.4	123.6	125.9	128.2	130.5	132.9	135.3	137.7	140.1	142.5
32	95.7	97.8	99.9	102.0	104.2	106.3	108.5	110.7	113.0	115.2	117.5	119.8	122.1	124.5	126.8	129.2	131.6	134.0	136.5	139.0	141.5	144.0	146.6	149.1
33	100.0	102.2	104.4	106.6	108.9	111.2	113.5	115.8	118.2	120.5	122.9	125.4	127.8	130.3	132.8	135.3	137.8	140.4	143.0	145.6	148.3	150.9	153.7	156.4
34	104.7	107.0	109.3	111.7	114.0	116.4	118.9	121.3	123.8	126.3	128.8	131.4	134.0	136.6	139.2	141.9	144.6	147.4	150.1	152.9	155.7	158.6	161.5	164.4
35	109.7	112.2	114.6	117.1	119.6	122.2	124.7	127.3	129.9	132.6	135.3	138.0	140.8	143.6	146.4	149.2	152.1	155.0	158.0	161.0	164.0	167.1	170.2	173.3
36	115.2	117.8	120.4	123.0	125.7	128.4	131.1	133.9	136.7	139.5	142.4	145.3	148.3	151.3	154.3	157.3	160.5	163.6	166.8	170.0	173.3	176.6	179.9	183.3
37	121.3	124.0	126.8	129.6	132.4	135.3	138.2	141.2	144.2	147.3	150.3	153.5	156.7	159.9	163.1	166.5	169.8	173.2	176.7	180.2	183.7	187.3	191.0	194.7
38	127.9	130.8	133.8	136.8	139.9	143.0	146.2	149.4	152.6	155.9	159.2	162.6	166.1	169.6	173.2	176.8	180.4	184.2	188.0	191.8	195.7	199.7	203.7	207.7
39	135.3	138.5	141.7	145.0	148.3	151.7	155.1	158.6	162.1	165.7	169.4	173.1	176.9	180.7	184.7	188.7	192.7	196.8	201.0	205.3	209.6	214.0	218.5	223.0
40	143.7	147.1	150.6	154.2	157.8	161.5	165.3	169.1	173.0	177.0	181.1	185.2	189.4	193.7	198.1	202.5	207.1	211.7	216.4	221.1	226.0	231.0	236.0	241.1
41	153.2	157.0	160.9	164.8	168.9	173.0	177.2	181.5	185.8	190.3	194.8	199.5	204.2	209.1	214.0	219.1	224.2	229.4	234.8	240.2	245.8	251.5	257.2	263.1
42	164.3	168.6	172.9	177.3	181.9	186.5	191.3	196.1	201.1	206.2	211.4	216.7	222.2	227.7	233.4	239.2	245.2	251.3	257.5	263.8	270.3	276.9	283.6	290.5
43	177.5	182.3	187.3	192.4	197.6	202.9	208.4	214.0	219.8	225.8	231.8	238.1	244.5	251.0	257.7	264.6	271.7	278.9	286.3	293.8	301.5	309.4	317.4	325.7
44	193.6	199.3	205.1	211.0	217.2	223.5	230.0	236.7	243.6	250.8	258.1	265.6	273.3	281.2	289.4	297.8	306.3	315.1	324.1	333.3	342.8	352.4	362.3	372.4
45	214.1	220.9	227.9	235.2	242.7	250.4	258.4	266.7	275.3	284.1	293.3	302.6	312.3	322.3	332.5	343.0	353.8	364.9	376.2	387.9	399.8	412.0	424.5	437.4
46	241.5	250.0	258.9	268.2	277.8	287.8	298.1	308.8	319.9	331.4	343.3	355.5	368.1	381.1	394.5	408.3	422.5	437.1	452.0	467.4	483.3	499.6	516.3	533.5
47	280.9	292.4	304.4	316.9	330.0	343.6	357.8	372.5	387.7	403.4	419.8	436.6	454.1	472.1	490.7	509.9	529.8	550.4	571.7	593.8	616.7	640.5	665.3	691.0
48	344.1	360.9	378.4	396.8	416.0	436.0	456.9	478.6	501.2	524.7	549.3	574.8	601.5	629.4	658.6	689.3	721.5	755.6	791.5	829.7	870.4	913.9	960.6	1011.2
49	461.1	488.4	517.2	547.5	579.4	613.1	648.8	686.7	727.0	770.1	816.4	866.4	920.8	980.4	1046.2	1119.9	1203.3	1297.7	1413.6	1553.1	1732.9	1966.3	2419.6	>2419.6

M

# Coliform Group

## Enzyme substrate detection reagent

### USER SPECIFICATIONS



**Character** This product is white or light yellow particles

**Clarification degree** Colorless or slightly yellow

**pH** 7.0 ~ 7.6

**Weight** 2.7 ± 0.5g

**Storage** Long term storage, drying, sealing and avoiding light storage at 4°C ~ 8°C

**Term of validity** 1 years

#### Working principle

In the water samples containing total coliform bacteria, the target bacteria were cultured in the ONPG-MUG medium at 36 ± 1 C. The specific enzyme betagalactosidase produced by the total coliform bacteria can decompose the color source substrate of the ONPG-MUG medium, which makes the culture medium yellow; meanwhile, Escherichia coli produces a specific beta-glucuronase to decompose the fluorescent substrate MUG in ONPG-MUG medium and produce characteristic fluorescence. The same principle, the heat tolerance coliform group (fecal coliform group) will decompose the color source substrate ONPG in the ONPG-MUG medium at 44.5 ± 0.5 °C, making the medium yellow.

**Packaging specification** 200 /box

#### Usage method

##### Qualitative detection

- Setp 1** Add reagent in 100ml water sample, dissolve it. and Culture 24h at 36 ± 1 °C
- Setp 2** Read result  
Colorless = negative  
Yellow = total coliform positive  
Yellow + fluorescence = Escherichia coli positive

Note: The heat tolerant coliform group (fecal coliform) needs to be cultured at 44.5 ± 0.5°C 24h, and yellow is positive.

##### Quantitative detection




- 1** Add reagent in 100ml water sample, dissolve it.
- 2** Pour into the 51 or 97 hole detection plate.
- 3** Use LK sealer to seal the detection plate, and cultivate 24h at 36 °C.
- 4** Results qualitative test was used to check the MPN table count.

Note: The heat tolerance coliform group (fecal coliform group) needs to be cultured at 44.5±0.5°C 24h, and the Yellow hole is positive, and the corresponding MPN table count is found.

Please read the instructions carefully before use

# Enzyme detection technology of Enterococcus

## USER SPECIFICATIONS



Please read the instructions carefully before use

**Character** This product is white or light yellow particles

**Clarity** Colorless or light yellow

**pH** 7.0 ~ 7.6

**Weight** 2.7 ± 0.5g

**Conditions** Storage at 4°C ~ 8°C, In cool dry place and protect from light.

**Validity** 1 Year, See reagent packaging for production date and expiration date.

### Science

Add water sample containing Enterococcus bacteria, culture the target bacteria in Mug medium at 41°C ± 0.5°C, and the specific biological enzymes produced by Enterococcus bacteria β-D-glucosidase can decompose the fluorescent substrate mug in the mug medium to produce β-D-glucoside (β-D-glucoside) and the characteristic fluorescent product 4-methyl-umbelliferone. Observe the fluorescence in the 366nm UV lamp, count through the quantitative detection disk, and query the MPN table to calculate the results.

**Package** 100 - test pack

### Usage method

#### Qualitative test

**Setp 1** Take 100ml water sample with 100ml sterile sampling bottle / quantitative bottle, add reagent, dissolve, and culture for 24h at 41°C ± 0.5°C.




**Setp 2** Result judgment  
 Non fluorescence = Negative results  
 Blue fluorescence = Positive results

#### Quantitative detection

- 1 Use 100ml sterile sampling bottle to measure 100ml water sample.
- 2 Add reagent, dissolve it, pour it into 51 Quanti-Tray or 97 Quanti-Tray.
- 3 Package the Quanti-Tray with a program-controlled quantitative sealing machine, and incubate it at 41°C ± 0.5°C for 24h.
- 4 Refer to the qualitative test for the result interpretation, and check the MPN table for counting.

# Enzyme detection technology of *Pseudomonas aeruginosa*

## USER SPECIFICATIONS



Please read the instructions carefully before use

**Character** This product is white or light yellow particles

**Clarity** Colorless or light yellow

**pH** 7.0 ~ 7.6

**Weight** 2.7 ± 0.5g

**Conditions** Storage at 4°C ~ 8°C, In cool dry place and protect from light.

**Validity** 1 Year, See reagent packaging for production date and expiration date.

### Science

Add water samples containing *Pseudomonas aeruginosa*, and pseudotech reagent contains a variety of essential amino acids, vitamins and other nutrients for *Pseudomonas aeruginosa*, providing excellent conditions for rapid growth and reproduction of *Pseudomonas aeruginosa* cells. The active *Pseudomonas aeruginosa* strain will secrete an enzyme, which can decompose the substrate in the reagent and produce substances that can emit blue fluorescence under ultraviolet light. Observe the fluorescence with a 366nm UV lamp, count with a quantitative detection disk, and query the MPN table to calculate the results.

**Package** 100 - test pack

### Usage method

#### Qualitative test

**Setp 1** Take 100ml water sample with 100ml sterile sampling bottle / quantitative bottle, add reagent, dissolve, and culture for 24h at 38°C ± 0.5°C.

**Setp 2** Result judgment  
 Non fluorescence = Negative results  
 Blue fluorescence = Positive results

#### Quantitative detection

**1** Use 100ml sterile sampling bottle to measure 100ml water sample.

**2** Add reagent, dissolve it, pour it into 51 Quanti-Tray or 97 Quanti-Tray.

**3** Package the Quanti-Tray with a program-controlled quantitative sealing machine, and incubate it at 38°C ± 0.5°C for 24h.

**4** Refer to the qualitative test for the result interpretation, and check the MPN table for counting.



# Multiple Enzyme Technology standard plate-count bacteria

## Scientific principles

The total bacterial count detection reagent uses enzyme substrate technology to detect the total bacterial count in water. The reagent contains a variety of unique enzyme substrates, each designed for different bacterial enzymes. When different enzyme substrates are decomposed by enzymes released by different bacteria, they release fluorescent groups. By observing the number of fluorescent cells under the ultraviolet lamp with the wavelength of 365 nm or 366 nm, the total value of colonies can be obtained by looking up the table.



## TEST STEP >>>



**1** Into a brown bottle containing solid reagent Add 1ml water sample and 9ml sterile water, Shake to dissolve completely.



**2** Pour the fully dissolved solution into the The center of the detection disk is detected.



**3** Turn the detection disk, All samples are automatically filled Check the hole grid of the disc.



**4** When all the holes are filled, Tilt detection panel, Adsorb the excess liquid into the sponge strip.



**5** Invert the detection disk, In  $36 \pm 1^\circ\text{C}$  incubator Culture for 48h.



**6** The number of fluorescent holes is digitally displayed under the UV lamp, The results are obtained by comparing with the MPN table.

Note: If you need to test 0.1ml, 0.01ml or smaller volume water samples, you can add water samples after dilution for testing.

## Colony count MPN key table (95% confidence interval)

Number of positive holes	Total number of colonies(MPN/ml)	Lower limit	Upper limit	Number of positive holes	Total number of colonies(MPN/ml)	Lower limit	Upper limit
0	<2	<3	<14	43	120	88	162
1	2	0.3	14	44	124	91	167
2	4	1	16	45	128	95	173
3	6	2	19	46	132	98	178
4	8	3	22	47	137	102	183
5	10	4	25	48	141	106	189
6	12	6	27	49	146	109	195
7	15	7	30	50	151	113	201
8	17	8	33	51	156	117	207
9	19	10	36	52	161	121	213
10	21	11	39	53	166	125	220
11	23	13	42	54	171	130	227
12	26	15	45	55	177	134	234
13	28	16	48	56	183	139	241
14	30	18	51	57	189	144	249
15	33	20	54	58	195	149	257
16	35	22	58	59	202	154	265
17	38	23	61	60	209	159	273
18	40	25	64	61	216	165	282
19	43	27	67	62	223	171	292
20	45	29	70	63	231	177	302
21	48	31	74	64	239	183	312
22	51	33	77	65	248	190	323
23	53	35	80	66	257	197	335
24	56	38	84	67	266	204	347
25	59	40	87	68	276	212	361
26	62	42	91	69	287	220	375
27	65	44	94	70	299	229	390
28	68	47	98	71	311	238	407
29	71	49	102	72	324	248	425
30	74	51	106	73	339	258	444
31	77	54	109	74	355	270	466
32	80	56	113	75	372	282	491
33	83	59	117	76	392	296	519
34	86	62	121	77	414	311	551
35	90	64	126	78	440	328	589
36	93	67	130	79	470	348	636
37	97	70	134	80	507	371	695
38	100	73	139	81	555	398	775
39	104	76	143	82	623	432	899
40	108	79	148	83	738	476	1146
41	112	82	152	84	>738	>476	>1146
42	116	85	157				

# PROGRAM-CONTROLLED AND QUANTITATIVE SEALER

## LK-2014A



➤ Darkroom (light box)

➤ Handheld UV Analyzer

➤ 100ml quantitative bottle

➤ 97-well quantitative  
detection plate  
51-well quantitative  
detection plate

### TEST STEP



#### STEP 1

Add the water sample to the 100ml sampling bottle, add the enzyme substrate detection reagent, tighten the cap and shake until completely dissolved.



#### STEP 2

Add the sample to the quantitative assay plate.



#### STEP 3

Use the program-controlled quantitative sealing machine to distribute and seal the quantitative detection disk.



#### STEP 4

Put it into a constant temperature incubator for cultivation. If the total coliform or Escherichia coli is detected, it should be cultivated at  $36\pm 1^{\circ}\text{C}$ . For detection of coliform bacteria, it should be cultured at  $44.5^{\circ}\text{C}$ .



#### STEP 5

Control the number of positive (yellow) grids on the colorimetric disc. If you need to detect Escherichia coli, you need to illuminate the yellow positive grids with a purple light to show blue fluorescence, which is Escherichia coli positive.



#### STEP 6

According to the type of the well plate of the experiment, check the MPN table of 51 wells or 97 wells respectively to obtain the number of coliform bacteria in 100 ml. Data are presented against MPN.

# PARAMETER

## Program-controlled quantitative sealing technical indicators

NO.	TECHNICAL INDEX	TECHNICAL PARAMETER
01	Name	Program-controlled quantitative sealing machine
02	Use	For the detection of total coliforms, Escherichia coli, fecal coliforms in water quality by enzyme substrate method
03	Reliability	No leaks, no holes
04	Stability	Can detect more than 40,000 samples, with a service life of more than 5 years
05	Convenience	On/off and reverse buttons, automatic stop function Digital display window, cleaning window
06	Fast	No need sterile room, 24h detection of total coliforms, Escherichia coli, fecal coliforms in water
07	Weight	≤16kg
08	Size	39*27*30cm
09	Preheat time	≤14 min
10	Noise	≤50dba
11	Housing temperature	≤40°C
12	Operating Voltage	AC 220V ±10%, 50Hz
13	Sealing speed	51 holes / 97 holes Quantitative detection disc sealing time 12 seconds/piece
14	Working temperature	-10°C~50°C
02	Detection range	With 51-hole quantitative detection plate detection range, water sample is not diluted. Compatible with the detection range of the 97-well quantitative detection plate, the water sample is not diluted.

Note: The equipment comes standard with 51-hole rubber pad and 97-hole rubber pad, operation manual, power cord, power fuse, and ISO and CE certification.

## Handheld Uv Analyzer With Darkroom

NO.	SPECIFICATIONS	PERFORMANCE
01	Name	Handheld UV Analyzer
02	Wavelength	Dual wavelength (254nm, 366nm)
03	Operating Voltage	AC 220V ±10%, 50Hz

## Consumables

NAME	MODEL	UNIT	QUANTITY	PACKAGING
Enzyme substrate method detection reagent	Colitech 24h	box	1 box	200 pcs/box
Quantitative bottle	with 100ml scale	box	1 box	200 pcs/box
Quantitative orifice plate/ Quantitative detection plate	51 or 97 holes	cartons	1 cartons	100 pcs/carton

# WATER TEST SOLUTION

NO.	ITEM NAME	SPECIFICATION	QUANTITY
<b>PROGRAM-CONTROLLED AND QUANTITATIVE SEALER</b>			
01	Host	/	1 pcs
02	Power cable	/	1 pcs
03	51-hole rubber gasket	/	1 pcs
04	97-hole rubber gasket	/	1 pcs
05	51 well MPN table	/	1 pcs
06	97-well MPN table	/	1 pcs
07	Certificate	/	1 pcs
08	Manual	/	1 pcs
09	Fuse	/	1 pcs
10	Sealer Tray	/	1 pcs
<b>51-well quantitative well plate/detection plate</b>			
11	Quantitative range 0-200MPN/100ml	100 pcs/carton	1 carton
<b>97-well quantitative well plate/detection plate</b>			
12	Quantitative range 0-2419MPN/100ml	100 pcs/carton	1 carton
<b>Colitech Enzyme Substrate Assay Reagents</b>			
13	Enzyme detection technology of Coliform group	200 pcs/box	1 box
	Enzyme detection technology of Enterococcus		
	Enzyme detection technology of Pseudomonas aeruginosa		
<b>100ml Sampling bottle with scale</b>			
14	Contains sodium thiosulfate	200 pcs/box	1 box
<b>UV Analyzer with Darkroom</b>			
15	366nm, 254nm	Dual Wavelength	1 pcs
<b>Color Pate</b>			
16	97-well positive cuvette		1 pcs
17	51-well positive cuvette		1 pcs
<b>Detection related tools (Provide by yourself)</b>			
Incubators, pipettes, tips and tip boxes, sterilizers, etc.			

DEVOTION TO LIFE IN COSMOS

## Lifecosc Biotech Limited

Add: 3B Fusen Building, 669 Shuangming Avenue, Guangming,  
Shenzhen, China

Tel:+86 755 23420845

Email: [info@lifecosc.com](mailto:info@lifecosc.com)

Web: [www.lifecosc.com](http://www.lifecosc.com)